



***1540 Mode**

1540 GRID RX (GFX)	12mj ↓ / 2 stack ↓	Downtime X
1540 Motion Toning (MT)	Moving	Downtime X
1540 DEEP FX (DFX)	12mj ↑ / 3 stack ↑	Downtime O

***1540 Treatment Guide**

Lesion	Mode	Energy	Distance	i-stack (Sub-pulse)	scanning	Pass	Tx size	Remark
Melasma	1540 GFX	8.4mj ~	0.9 ~ 1.0mm	1 ~ 2nd	Random	1 pass	□ 15 x 15	Toning treatment program +α
	1540 MT	8.4mj	0.6mm	1st	i-motion	2 pass	○ 15 x 15	
Rejuvenation (Skin tone/ texture & Mild pore)	1540 GFX & 1540 MT	10.8 ~ 12.0mj	0.9mm	2nd	Random	1 pass	□ 15 x 15	
		10.8 ~ 12.0mj	0.6mm	1st	i-motion	5 pass	○ 15 x 15	
	Only 1540 MT	10.8 ~ 12.0mj	0.6mm	1st	i-motion	10 pass	○ 15 x 15	
	1540 DFX	12.0mj ~	0.9mm	3rd	Random	1 pass	□ 15 x 15	
Hair Regrows	1540 GFX	9.6 ~ 10.8mj	0.8mm	2nd	Random	1 ~ 2 pass	□ 15 x 15 ⇐ 15 x 15	10 times every week + Solution, + Secret
Acne	1540 GFX (P)	10.8mj ~	0.6mm	2nd	Random	1 ~ 2 pass	○ 6 x 6	Point treatment for the purulent area
	1540 MT	10.8mj ~	0.6mm	1st	i-motion	4 ~ 5 pass	○ 15 x 15	
Acne Scar & Scar	1540 DFX	13.2mj ~	0.6mm	3rd	Random	1 pass	Set it according to the treatment area	Point treatment for the treatment area

***1540nm FX Laser Single Treatment (1540 GFX & 1540 MT)**

- A program that leads patients through repetitive procedures without using consumables and without downtime.

***1540 Single Program**

Lesion	Treatment Intervals and Programs					
Melasma	Treatment is performed by including as +α in the toning program currently being used in the clinic. (1540 GFX & 1540 MT) This procedure is Fractional Thermal Toning, which induces dermal regeneration by delivering appropriate heat to the target area while minimizing damage to the epidermis.					
Rejuvenation (Skin tone/ texture & Mild pore)	① 1540 GFX + 1540 MT + Skin booster Solution ② 1540 MT + Skin booster Solution		The procedure is performed 5 times as a standard, and you can decide whether to add 5 times according to the skin condition of the patient.			
Hair Regrows	① 1540 GFX or 1540 MT + Scalp Solution		The procedure is performed five times every week, and the procedure is added five times depending on the patient's skin condition.			
Acne	1st procedure	2nd procedure	3rd procedure	4th procedure	5th procedure	+α
	1540 GFX(Point treatment) & MT	1540 MT	1540 GFX & MT	1540 MT	1540 MT	In case of necessary
	The treatment interval is once a week, and after 5 treatments, you can decide whether to proceed with additional treatment according to the patient's condition.					