



Secret DUO Combination Treatment / Scar / 1 session

Mode	Pass	Guideline		
RF Micro-needling	1 st Pass ~2 nd Pass	<p>Rolling scar – 25 Pin Insulated needle, - Use it to burn tissues which holds the skin to be pulled down wards - To coagulation at dermis (deeper layer)</p> <p>1st 60~70% / 100ms / 2mm 2nd 60~70% / 100ms / 1.5mm</p> <p>Box scar, ice peak scar - 25 Pin Non-insulated needle epidermal stimulate, do not use high parameter, it might be burned the skin</p> <p>1st 50% / 150ms / 2mm 2nd 40% / 200ms / 1.5mm</p>		
1540 Erbium Glass Laser	3 rd Pass ~4 th Pass	<p>3rd 300um lens Energy 15.6mj / i-Stack 1st (Or 8.4mj/i-Stack 2nd)</p> <p>Distance 0.8mm Scanning size 10x10 Shape : Circle Scanning : I-MOTION Overlap : 50%</p> <p>4th Same as 3rd or change energy level up to 18mj / i-Stack 1st</p>	<p>Skip keep jumping (ex.jump every 2 beeps sound)</p> <p>: To keep enough time while the laser penetrates fully since pulse duration is 1300us (based on 15.6mj Energy)</p> <p>: Fast moving not be able to finish destored before we move.</p>	<p>i- Motion Scanning : Can keep the energy density evenly at the overlap space</p> <p>i- Motion</p>  <p>Random</p> 

- Based on Asian skin type. please proceed pre-test before treatment
- Provided PT. Zhengan Indonesia

Secret DUO 1540nm Laser Treatment / Scar, Acne, Skin Toning

Indication	Applicator	Protocol	Remark
SCAR	300um Lens	Energy 8.4mj (up to 10.8mj) Distance 0.8mm I-stack 3 rd Scanning size : 10x10 Shape : Circle Scan type : i-Motion 2 passes	Give the higher damage to make a coagulation columns into the skin, Stamping movement. Downtime required If required more higher energy value, can increase the energy up to 10.8mj, do not over I-stack 3 rd .
ACNE	600um Lens	Energy 10.8 ~ 12mj Distance 0.6mm I-stack 1st Scanning size : 10x10 Shape : Circle Scan type : i-Motion 2~4 passes	Quickly move as like Q-switched Nd:YAG Toning, Effect of the heat accumulation
SKIN TONING	600um Lens	Energy 13.2 ~ 15.6mj Distance 1.0mm I-stack 1st Scanning size : 15x15 Scan type : i-Motion 2 passes –up to 4 passes Only need to slight overlap	- No anesthesia, Quick treatment and Immediate effect, painless treatment, NO Downtime - For Brighter Skin, smooth texture lasting 1 month up to 2months

- Based on Asian skin type. please proceed pre-test before treatment
- Provided PT. Zhengnan Indonesia